AQUALARM®

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CE 0459

What is Aqualarm and when is it used?

Aqualarm is a carbomer- and ion-based fluid ophthalmic gel which lubricates the eye surface, i.e. it creates a transparent, lubricating, moisturizing film on the eye surface, providing temporary compensation for a lack of tears.

Carbomer forms an even, stable, transparent gel with good "at rest" surface coverage. It spreads easily across the eye surface, creating an effective barrier for the epithelium (cells on the surface of the eye). It stays on the eye surface for a long time. Longer contact time means less frequent need for treatment.

pH and osmolality are similar to those found in tears.

Aqualarm is used to give relief to dry, tired eyes, and can also be used to combat other minor eye problems connected with smoke, pollution, air-conditioning, dust, dry heat or prolonged work in front of a screen.

What does Aqualarm contain?

Carbomer 980: 2 mg per 1g of solution.

Povidone K30, disodium edetate, sodium hydroxide, monosodium phosphate dihydrate, disodium phosphate dodecahydrate, sorbitol, purified water.

The solution is preservative-free. It is sterilised by autoclaving.

How is Aqualarm used?

Instructions for use:

- Detach the single-dose container and open it by twisting off the removable tab.
- Look upwards, pull the lower eyelid slightly downwards and place one drop of solution in the inner corner of the eye. Repeat 4-5 times a day.

Precautions:

- Do not touch the eye with the single-dose container nozzle.
- Wash hands thoroughly before use.
- If Aqualarm is used in association with other ophthalmic products, wait at least 5 minutes between the two treatments and always use Aqualarm last.

When should Aqualarm not be used?

If there is allergy to one of its components.

What precautions should be taken when using Aqualarm?

Aqualarm is a single-dose solution, the single-dose container must be thrown away immediately after use.

Do not use a single-dose container which has already been opened and/or partially used.

In case of re-use, the laboratory does not guarantee cleanness and sterility of the product.

Do not use after the expiry date.

Do not use if the container is damaged.

The product should not be injected or taken by mouth.

If vision is impaired, you are advised not to drive or use machinery.

Keep out of the reach of children.

How should the product be stored?

Store at a temperature below 30°C

Laboratoire CHAUVIN
416, rue Samuel Morse
34000 MONTPELLIER - France

AQUALARM is a registered trademark of
Bausch & Lomb Incorporated or its affiliates
©Bausch & Lomb Incorporated

CE 0459